

MANUAL HANDLING

Walk Transfer Belts



Walk Transfer Belts

Exercise is an important part of the patient's therapy, whether it be part of a rehabilitation programme, or keeping the independent client mobile. The Walk Transfer Belt was designed to promote walking as therapy, whilst providing comfort and security. When the Walk Transfer Belt is worn around the waist, the client is able to be stabilised and guided by the carer using the hand holds on the belt.

The Walk Transfer Belt can also be used as a transfer belt for stand transfer or positional guidance. Walk/Transfer belts have optional rubber backing for extra grip to client and are made with padded nylon.

The Walk Transfer Belt is tested to International Standards at 150kg.



Cam Lock Buckle



1 Way Side Release Buckle

Features:

- Soft, strong nylon padded belt for extra comfort
- PVC padded belt for ease of cleaning between clients
- Optional rubber backing on nylon material belts to prevent slipping
- Adjustable for snug fit
- Various closures to choose from: Velcro, Side Release Buckle, Cam Lock Buckle
- Four padded hand holds allow variation for position of carers
- Fully washable
- Used to assist client from sitting to standing position (or vice-versa)
- Used to stabilise and support a client whilst walking, either with or without a walking frame (depending upon client's condition)

Sizing:

- X-Small fits waist 50cm - 75cm
- Small fits waist 65cm - 95cm
- Medium fits waist 80 cm - 110cm
- Large fits waist 100cm - 130cm
- X-Large fits waist 125cm - 155cm

Types of Belts:

- MHTA-WTB1B: 1 Way Side Release Buckle Closure
- MHTA-WTB1BR: 1 Way Side Release Buckle Closure /Rubber Backing
- MHTA-WTB2B: 2 Way Side Release Buckle Closure
- MHTA-WTB2BR: 2 Way Side Release Buckle Closure /Rubber Backing
- MHTA-WTBCL: Cam Lock Buckle Closure
- MHTA-WTBCLR: Cam Lock Buckle Closure / Rubber Backing
- MHTA-WTBV: Velcro Closure
- MHTA-WTBVR: Velcro Closure / Rubber Backing

SHA221014-3

MANUAL HANDLING

Walk Transfer Belts



Product Code	Product Name	Size
MHTA-WTB1B-XS	Walk Transfer Belt: 1 x Way Side Release Buckle Closure	X-Small
MHTA-WTB1B-SM	Walk Transfer Belt: 1 x Way Side Release Buckle Closure	Small
MHTA-WTB1B-MD	Walk Transfer Belt: 1 x Way Side Release Buckle Closure	Medium
MHTA-WTB1B-LG	Walk Transfer Belt: 1 x Way Side Release Buckle Closure	Large
MHTA-WTB1B-XL	Walk Transfer Belt: 1 x Way Side Release Buckle Closure	X-Large
MHTA-WTB1BR-XS	Walk Transfer Belt: 1 x Way Side Release Buckle Closure/Rubber Backing	X-Small
MHTA-WTB1BR-SM	Walk Transfer Belt: 1 x Way Side Release Buckle Closure/Rubber Backing	Small
MHTA-WTB1BR-MD	Walk Transfer Belt: 1 x Way Side Release Buckle Closure/Rubber Backing	Medium
MHTA-WTB1BR-LG	Walk Transfer Belt: 1 x Way Side Release Buckle Closure/Rubber Backing	Large
MHTA-WTB1BR-XL	Walk Transfer Belt: 1 x Way Side Release Buckle Closure/Rubber Backing	X-Large
MHTA-WTB2B-XS	Walk Transfer Belt: 2 x Way Side Release Buckle Closure	X-Small
MHTA-WTB2B-SM	Walk Transfer Belt: 2 x Way Side Release Buckle Closure	Small
MHTA-WTB2B-MD	Walk Transfer Belt: 2 x Way Side Release Buckle Closure	Medium
MHTA-WTB2B-LG	Walk Transfer Belt: 2 x Way Side Release Buckle Closure	Large
MHTA-WTB2B-XL	Walk Transfer Belt: 2 x Way Side Release Buckle Closure	X-Large
MHTA-WTBCL-XS	Walk Transfer Belt: Cam Lock Buckle Closure	X-Small
MHTA-WTBCL-SM	Walk Transfer Belt: Cam Lock Buckle Closure	Small
MHTA-WTBCL-MD	Walk Transfer Belt: Cam Lock Buckle Closure	Medium
MHTA-WTBCL-LG	Walk Transfer Belt: Cam Lock Buckle Closure	Large
MHTA-WTBCL-XL	Walk Transfer Belt: Cam Lock Buckle Closure	X-Large
MHTA-WTBCLR-XS	Walk Transfer Belt: Cam Lock Release Buckle Closure/Rubber Backing	X-Small
MHTA-WTBCLR-SM	Walk Transfer Belt: Cam Lock Release Buckle Closure/Rubber Backing	Small
MHTA-WTBCLR-MD	Walk Transfer Belt: Cam Lock Release Buckle Closure/Rubber Backing	Medium
MHTA-WTBCLR-LG	Walk Transfer Belt: Cam Lock Release Buckle Closure/Rubber Backing	Large
MHTA-WTBCLR-XL	Walk Transfer Belt: Cam Lock Release Buckle Closure/Rubber Backing	X-Large
MHTA-WTBV-XS	Walk Transfer Belt: Velcro Release Buckle Closure	X-Small
MHTA-WTBV-SM	Walk Transfer Belt: Velcro Release Buckle Closure	Small
MHTA-WTBV-MD	Walk Transfer Belt: Velcro Release Buckle Closure	Medium
MHTA-WTBV-LG	Walk Transfer Belt: Velcro Release Buckle Closure	Large
MHTA-WTBV-XL	Walk Transfer Belt: Velcro Release Buckle Closure	X-Large
MHTA-WTBVR-XS	Walk Transfer Belt: Velcro Release Buckle Closure/Rubber Backing	X-Small
MHTA-WTBVR-SM	Walk Transfer Belt: Velcro Release Buckle Closure/Rubber Backing	Small
MHTA-WTBVR-MD	Walk Transfer Belt: Velcro Release Buckle Closure/Rubber Backing	Medium
MHTA-WTBVR-LG	Walk Transfer Belt: Velcro Release Buckle Closure/Rubber Backing	Large
MHTA-WTBVR-XL	Walk Transfer Belt: Velcro Release Buckle Closure/Rubber Backing	X-Large