

WHS

StaminaLift TS5000 Transfer System



Proven to prevent injuries related to moving hospital beds and equipment

The result of extensive research, development and field testing, our **StaminaLift TS5000 Transfer System** is unlike any other hospital bed mover. The TS5000 is easy to operate, and you will appreciate it for its universal compatibility – the ability to connect to a wide range of beds, stretchers and equipment.

**On
NSW Govt
Contract
217A**



With the TS5000 Transfer System you can move:

- A wide range of acute care hospital bed and stretchers including **Hill-Rom, ArjoHuntleigh, Stryker, Linet, Howard Wright, Prius, Volker** and more
- Some patient care chairs, anaesthetic machines, linen carts, trolleys, rubbish skips and more

Standards Compliance & Testing:

Complies with the relevant requirements of the Applicable Standard EN60601-1 (3rd Ed)

SHA110814-6

WHS

StaminaLift TS5000 Transfer System



Features:

- **Universal Jaw Assembly:** Designed for global acute care bed compatibility. Accommodates many stretchers and specialty beds
- **Adjustable Jaws:** Wheel base width: 550 - 890mm Castor sizes: 125mm - 200mm (single & double)
- **Adjustable Handle:** Comfort regardless of operator height
- **Sealed Gel Batteries:** 15 hours normal use before recharge
- **Variable Speed:** Flexibility for operator control
- **Battery Charger:** 120 or 240V for the global market

Success Story:

“In 2009, to address a high rate of wardsmen injuries, we introduced StaminaLift bed movers at our 855 bed hospital. There was an immediate and sustained decrease in injury rate and severity.”

- Leslie Gibbs, Work Health and Safety Manager
Liverpool Hospital, NSW Australia

Technical Specifications

Max load capacity	300 kg lift and 600 kg push or pull (370kg / 740 kg option on request)
Max incline capacity	7 degrees
Weight	160 kg
Dimensions	In use: W 780 mm x L 925 mm x H 1220 mm
Batteries	Sealed maintenance free GEL batteries. Up to 15 hours operation from full charge
Recharge	<ul style="list-style-type: none"> • Nightly recharge recommended using supplied automatic charger • Up to 6 hours to fully charge
Warning device	Horn (warning device) on control unit
Security	Electronic locking function available, refer to operators manual
Wheels	Non-marking, suitable for indoor use on a variety of surfaces
Controls	Programmable, ergonomic joystick for control of direction, speed and jaw movement
Safety features	<ul style="list-style-type: none"> • Automatic brakes / emergency shut-down button • Emergency jaw release – allows override of lowering mechanism • Emergency Brake release – allows override of brake mechanism • Safety features to prevent unauthorized use
Motor	Two variable drive DC electric motors

Optional StaminaLift Ride-On Platform Now Available!



StaminaLift Ride-On Platform
folded down



StaminaLift Ride-On Platform
folded up



StaminaLift Ride-On Platform
in use

Product Code	Product Name
OHSL-TS5000	StaminaLift TS5000 Transfer System
OHSL-S00105	StaminaLift Ride-On Platform